ARCEOLE

CE 0459

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Set for infusion of ophthalmic gas

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| INSTRUCTIONS FOR USE | | | | |
|--|----------------|------------------------------------|------------------------------------|---|
| ARCEOLE • INDICATIONS : Endotamponador de retina. Endotamponado de la córnea in Descemet membrane detachment and lamellar keratoplasty. | | | | |
| DESCRIPTION OF THE MEDICAL DEVICE | | | | |
| ARCEOLE is composed of a sterile single-dose cylinder of 30ml non sterile gas for ophthalmic surgery and a sterile plastic graduated 50ml syringe mounted with a sterilizing 0.22µm filter and a device to connect it onto the cylinder. The measuring tolerance of the 50 ml syringe is 4% of the measured volume. A 30G needle is supplied for the intraocular injection. A patient wrist band allows the identification of intraocular gas wearer. A set of 5 patient labels is provided to allow traceability. SET PS and SET 1 ml are sold separately. | | | | |
| The different accessories provided with ARCEOLE have been validated for the hereunder detailed use. It is not advisable to substitute them by others. | | | | |
| Nota : SET PS and SET 1 ml se presentan en blister estériles containing 1 syringe mounted with a 0.22µm filter and a connector (SET PS: 5 ml syringe, SET 1ml: 1 ml syringe). The must be used according to the directions indicated below. The measuring tolerance of the 5 ml syringe is 4% of the measured volume and 5% of the measured volume for the 1 ml syringe. | | | | |
| COMPOSITION | | | | |
| ARCEOLE is a sterile single-dose cylinder of one of the following gases: SF6 (Sulfur hexafluoride), C2F6 (Hexafluoroethane), C3F8 (Octafluoropropane). | | | | |
| Each sterile single-dose cylinder contains a pure non sterile gas in a sterile packaging. Does not contain latex or phthalates | | | | |
| Nota : The gases included in ARCEOLE are non-toxic, inert, non-flammable, colorless and odorless compounds. | | | | |
| Once implanted, they are not metabolized and gradually eliminated through bloodstream and the respiratory tract. | | | | |
| Compatibilité with MRI : ARCEOLE is composed of perfluorinated compounds which do not contain hydrogen, they therefore generate no signal. | | | | |
| Expansion and re-absorption rates based on bibliography data: | | | | |
| Gas | Expansion rate | Time allowed for expansion in days | Average tamponade duration in days | Non expansive concentration in % |
| SF6 | 1,9-2 | 2 | 10-15 | 20 |
| C2F6 | 3,3 | 3 | 30-35 | 16-17 |
| C3F8 | 4 | 3 | 55-65 | 12 |

CONTRAINDICATIONS

It is imperative to stop the anaesthesia with nitrogen protoxide inhalation at least 15mn before starting the gas injection into the patient eye.

A patient with a gas bubble must absolutely avoid nitrogen protoxide anaesthesia during the gas bubble remaining period and at least for 3 months following the date of implantation. Severe post-operative complications up to blindness linked to the intraocular pressure increase may occur.

In order to draw attention to those precautions, the patient needs to wear the wristband included in the packaging.

A patient implanted with intraocular gas must absolutely avoid pressure variation during the gas bubble remaining period and at least for 3 months following the date of implantation (elevation, travel by plane, diving with or without isobar caisson, hyperbaric caisson treatment). The resulting intraocular pressure increase may lead to blindness.

PRECAUTIONS AND WARNINGS

The gases for ophthalmic surgery are exclusively intended to be used by ophthalmologists trained to this technique.

Although not toxic these gases are competing with oxygen when inhaled. Avoid intentional inhalation.

A preoperative gonioscopy is recommended! Avoid injection of pure gas in case of narrow irido-corneal angle.

Before and after injection of gas, monitoring of the retina central artery is required.

If an expansive gas is used, it is recommended to monitor the intraocular pressure during the gas expansion phase.

It is recommended to proceed with a suture at the end of the 23 G transconjunctival vitrectomy to maintain the eye gas-tight. In case of a large Descemet membrane tear, slowly inject the gas with the canula oriented in the iris plane to avoid fish egg formation.

After tamponade, a treatment against ocular hypertension should be established and IOP daily monitored.

Perform regular follow up of glaucoma patients or patients with IOP, the elderly patients and those with diabetes, patients with inferior tear which may increase a secondary tear risk, peroperative haemorrhage risk, and patients with indentation.

As a rule, it is necessary to check the ability of the patient to hold his head in the required position during the treatment.

The patient head positioning requirement may be waived whrm the posterior chamber is filled up to 70% in the case of vitrectomy. The risk of appearance of cataract latérale is to be avoided.

In **vitreo-retinal surgery**, **ARCEOLE** is a medical device designed to inject an ophthalmic gas mixed with air. In case of injection of pure gas, under the physician responsibility, the latter shall take into account the increased volume of the bubble of the pure gas compared to a non-expansive gas bubble. Toxicity in the Rabbit has been shown in relation with the dose and tamponade duration.

COMPLICATIONS AND SIDE EFFECTS

The complications described in the literature are those linked to ophthalmic surgery and particularly:

Vitreo-retinal surgery

- Sub-retinal passage, hemorrhage, haematoma, choroid tear, sub-conjunctival passage, anterior chamber passage (aphakic patients)

- Risky lens or retina touch with the needle during the injection.

- Risk of separation or capture of the gas bubble between the crystalline lens and the pars plana.

- Refractive change

- Anterior chamber flattening

- Early gas bubble decrease in the case of self-sealing incision

- Extension of the tear to the macula, inadequate position of the gas bubble (respect the patient head positioning during the injection)

- Vitreous change

- Retinal or vitreal capture at te injection site

- Fogging on intraocular lenses especially silicone lenses.

Potential side effects are: central artery occlusion, hypertonia, secondary retinal tear, cataract, inflammation, endophtalmis, keratopathy, hemorrhage, vitreous reaction, macula hole and must be considered by the surgeon with regards to the expected benefits.

Corneal surgery

Hypertonia, risk of gas bubble fragmentation or capture between the parts to reattach, infection, pupillary block.

In the case of irita-corneal injection, caution is advised when the patient is implanted with an hydrophilic acrylic intraocular lens, a case of unexplained opacification has been reported in the bibliography (*SAEED MJ, SINGH AJ, MOREL AJ, Sequential Descemet's membrane detachments and intraocular lens haze secondary to SF6 or C3F8. Eur. J. Ophthalmol. 2006; 16: 758-60*).

INSTUCTIONS FOR USE

WARNING: the gas is not sterile and must be sterilized by passing through the sterilizing filter supplied with the product. Preparation must be performed extemporaneously.

Before using the syringe, it is advised to release the plunger stopper from the bottom of the syringe barrel by pulling the plunger. Remove the needle from the barrel and smoothly push the plunger to express the air from the syringe barrel.

Remove the filter from the syringe.

Immediately mount the sterile 30G injection needle with its protective cap, the mix is ready for intraocular injection.

Pure gas

As an indication and according to the bibliography, the volume to be injected is between 0.3 to 0.6 ml for the SF6, 0,4 ml for the C2F6, and 0,3 ml for C3F8 (*MATHIAS A, PAGOT V. Tamponnements internes EMAC (Paris-France) Ophthalmologie, 1995 ; 21 248A –29, 10p*).

Corneal surgery

In the lamellar keratoplasty and Descemet membrane reattachment, the amount of gas to be injected is about 0,25 ml, preferably mixed with 20% air for SF6 and 12% air for C3F8.

SET 1 ml is recommended for injection of the ophthalmic gases in the anterior chamber.

The table below provides help to define the gas concentration:

| Gas | Non expansive concentration in % | Volume of pure gas to fill a 1ml syringe | Volume of air |
|------|---|--|---------------|
| SF6 | 20 | 0,20 ml | 0,80 ml |
| C3F8 | 12 | 0,12 ml | 0,88 ml |

Purge the SET 1 ml syringe from its gas content one time and re-fill as indicated above.

Take air in the syringe to obtain the desired proportion mix. The air is sterilized by passing through the syringe filter.

Remove the filter from the syringe.

Immediately mount the 30G sterile injection needle with its protective cap, the mix is ready for intraocular injection.

STORAGE CONDITIONS

Store in dry, aerated area and at temperature <25°C.

Do not use if the sterile pouch has been compromised.

Do not use after expiration date.

Single use. Purity and sterility of the product are guaranteed for a single use only.

Do not re-sterilize.

WASTE DISPOSAL:

| With household waste or recycle | With hospital waste contaminated products |
|---|---|
| Paper Box, instructions for use leaflet and traceability labels | |
| PVC hard blister, HDPE blister lid | 30G needle filter |
| Aluminium canister (no practical residual pressure, no flammable gas) | |
| 50 ml syringe, 5 ml syringe, 1 ml syringe, connector : polypropylene | |

STERILIZATION

The gas contained in **ARCEOLE** is not sterile, only the canister and accessories are sterilized with ethylene oxide.

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