



7900

Ophthalmic Instrument Stand

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



7900 Ophthalmic Instrument Stand

Contemporary Design Matches Reliance Chairs

The 7900 IC (with instrument charging console) and the 7900 NC (no console) provide maximum access to patients and are compatible with all Reliance chairs. They include chair controls and three wells for handheld instruments (IC version).

Arms for Instrument Delivery

Counterbalanced slit-lamp arm locks and allows adjustment with a single hand. 12" of vertical travel can be counterbalanced for instruments 12 to 60 lbs. A raise and lower switch on the slit lamp arm controls chair height adjustments.

Upper refractor arm offers 10" range of vertical movement; delivers instruments up to 20 lbs.

Third instrument arm available as an option: 10- to 50 lb. capacity; 10" of vertical travel.

5399 Monitor Arm: Instrument stand supports optional 5399 monitor arm.

Durability Features

Space-age base material including the same polycarbonate used on supersonic aircraft, resists scrapes, scuffs and scratches.

Backlit membrane switches last 10-times longer than mechanical switches and increase safety with lower voltages. They control the chart projector, overhead lamp and binding post for an indirect ophthalmoscope, and raise and lower the exam chair.

Halogen lamp with long-lasting bulb provides bright, direct light for minor procedures. Arm adjusts without tools. Low voltage for safety.

CSA mark of quality indicates that the 7900 has passed all rigorous CSA, C/US tests, affirming its safety and durability and achieving certificate of compliance.

OPTIONS

Wheelchair Access

Deliver care to wheelchair bound patients with an optional lower slit lamp arm. Cost applies to Americans With Disability Act's Tax Credit.

Eclipse Room Light Control

The Eclipse activates all room lighting from the stand, or a small console. Lighting adjusts based on instrument moves or quick manual adjustments.

Tech Specs

Height to console: 36". Height to pole top: 71". Base: 21x 21". Shipping weight: 292 lbs. Electrical (Standard UL Listed), 115 volts, 60 Hz, 3 amperes. Optional: 230 volts, 60 Hz, 1.5 amperes.

