





Anaesthetic Cannulas and Needles



















Anaesthetic Cannulas

For administering anaesthetic into the sub tenon space.

M0278	19 G		19 G x 1" 1.1 x 25mm	Sub tenon cannula. Curved	Flattened tip		BOX	
							10	100
M0279	20 G		20 G x 1" 0.9 x 25mm	Sub tenon cannula. Curved	Flattened tip		10	100

Anaesthetic Needles

For administering anaesthetic around the globe.

M0275	25 G		25 G x 1½" 0.5 x 38mm	Retrobulbar. Tri-facet		BOX	
						10	100
M0637	25 G		25 G x 1½" 0.5 x 38mm	Retrobulbar (Atkinson)		10	100
M0638	23 G		23 G x 1½" 0.64 x 38mm	Retrobulbar (Atkinson)		10	100
M0641A	27 G		27 G x 1¼" 0.4 x 32mm	Peribulbar (Atkinson)		10	100
M0642	25 G		25 G x 1¼" 0.5 x 32mm	Peribulbar (Atkinson)		10	100
M0643	23 G		23 G x 1¼" 0.64 x 32mm	Peribulbar (Atkinson)		10	100
M0641	27 G		27 G x ⅞" 0.4 x 22mm	Peribulbar (Atkinson)		10	100
M0642A	25 G		25 G x ⅞" 0.5 x 22mm	Peribulbar (Atkinson)		10	100