




Ophthalmic Drape Portfolio

BVI offers a wide variety of options to meet your clinical needs



Half Body Drapes

Plastic Incise Drape with Fluid Collection Bag

Qty: 10/Bx

Description		SKU #
Plastic Incise Drape with Fluid Collection Bag		581149
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Polyethylene Opaque Film	



Plastic Aperture Drape with Fluid Collection Bag

Qty: 10/Bx


Description		SKU #
Plastic Aperture Drape with Fluid Collection Bag		581151
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Polyethylene Opaque Film	



Non-Woven Incise Drape with Fluid Collection Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Incise Drape with Fluid Collection Bag		581150
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	




Half Body Drapes

Non-Woven Aperture Drape with Fluid Collection Bag

Qty: 10/Bx


Description		SKU #
Non-Woven Aperture Drape with Fluid Collection Bag		581152
Overall size:	100 cm x 127 cm (39 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag:	18 cm x 16.5 cm (7 in x 6.5 in)	
Capacity:	250 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Incise Drape with Channel Bag

Qty: 10/Bx


Description		SKU #
Non-Woven Incise Drape with Channel Bag		581153
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	



Non-Woven Aperture Drape with Channel Bag

Qty: 10/Bx

Description		SKU #
Non-Woven Aperture Drape with Channel Bag		581155
Overall size:	140 cm x 127 cm (55 in x 50 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Fluid collection bag with channels:	26 cm x 35 cm (10 in x 14 in)	
Capacity:	500 ml	
Material Type:	Non-Woven Wetlaid	




Half Body Drapes

Non-Woven Half Body Incise Drape* with Channel Bag

Qty: 10/Bx

Description	SKU #
Non-Woven Half Body Incise Drape* with Channel Bag	581166
Overall size:	160 cm x 180 cm (63 in x 71 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Channel bag:	62 cm x 40 cm (24 in x 16 in)
Capacity:	1000 ml
Material Type:	SMS

*Without wick




Full Body Drapes

Full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Full Body Incise Drape

Qty: 10/Bx

Description	SKU #
Non-Woven Incise Drape with 2 Bags	581156
Overall size:	140 cm x 244 cm (55 in x 96 in)
Incise area:	13 cm x 15 cm (5 in x 6 in)
Fluid collection bag:	25 cm x 15 cm (10 in x 6 in)
Capacity:	2 x 500 mL
Material Type:	Non-Woven Wetlaid




Full Body Drapes

Full Body Aperture Drape

Qty: 10/Bx


Description		SKU #
Full Body Aperture Drape with 2 Bags		581157
Overall size:	140 cm x 244 cm (55 in x 96 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Capacity:	2 x 500 mL	
Material Type:	Non-Woven Wetlaid	



Full Body Incise Drape with Plastic Head and Channel Bag

Qty: 10/Bx

Description		SKU #
Full Body Incise Drape with Plastic Head and Channel Bag		581163
Overall size:	160 cm x 252 cm (63 in x 99 in)	
Incise Area:	13 cm x 15 cm (5 in x 6 in)	
Channel Bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	SMS 157 cm (61 in), Plastic 100 cm (39 in)	




*Without wick

Non-Woven Full Body Incise Drape

Qty: 10/Bx

Description		SKU #
Non-Woven Full Body Incise Drape		581164
Overall size:	160 cm x 252 cm (60 in x 99 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Channel bag:	62 cm x 40 cm (24 in x 16 in)	
Capacity:	1000 ml	
Material Type:	Non-Woven Wetlaid	



*Without wick

Face Drapes

Ideal for refractive or office procedures.


Mini Pre-Split Incise Drape

Adhesive backing.

Qty: 10/Bx

Description		SKU #
Mini Pre-Split Incise Drape		581154
Overall size:	40 cm x 40 cm (16 in x 16 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Pre-split:	5 cm (2 in)	
Material Type:	SMS	

*Without wick



Mini Plastic Incise Drape

Adhesive backing.

Qty: 10/Bx

Description		SKU #
Mini Plastic Incise Drape		581159
Overall size:	40 cm x 40 cm (16 in x 16 in)	
Incise area:	13 cm x 15 cm (5 in x 6 in)	
Material Type:	Polyethylene Opaque Film	

*Without wick



Mini Plastic Aperture Drape

Adhesive backing.

Qty: 10/Bx

Description		SKU #
Mini Plastic Aperture Drape		581160
Overall size:	40 cm x 40 cm (16 in x 16 in)	
Aperture area:	4 cm x 6 cm (1.5 in x 2 in)	
Material Type:	Polyethylene Opaque Film	

*Without wick



Dual Refractive Drape

Qty: 10/Bx

Description		SKU #
Dual Refractive Drape		581162
Overall size:	40 cm x 26 cm (16 in x 10 in)	
Aperture area:	2.5 cm x 4 cm (1 in x 1.5 in)	
Material Type:	Polyethylene Opaque Film	

*Without wick



Drape Accessories

Visitec® Eye Fluid Wick

Wicks absorb and direct fluid from surgical area into fluid collection bag.




Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Eye fluid wick, 15 cm Long x 0.4 cm Wide	581147

Mini Plastic Aperture Drape

Small incise provides complete coverage of lashes, lids and brows. Converts aperture to incise drape.





Qty: 10/Bx

Description	SKU #
Incise Drape	581161
Overall size:	10cm x 20 cm (4 in x 8 in)
Adhesive area:	11 cm x 16 cm (4 in x 6 in)
Material Type:	Polyethylene Opaque Film

Visitec® Fluid Collection Bags

For procedures where collection and control of fluid is designed to keep surgeon, patient and floor dry. Reduces operating room cleanup and turnaround time. Bags have an adhesive backing to adhere to surgical drape.



Qty: 20/Bx
Brand: Visitec®

Description	SKU #
Fluid Collection Bag: 17 cm x 16 cm (6.75 in x 6 1/4 in), 250 ml capacity	581177
Fluid Collection Bag with Wick Bag: 17 cm x 16 cm (6.75 in x 6.25 in), 250 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581146
Mini Fluid Collection with Wick Bag: 9 cm x 16 cm (3.5 in x 6.25 in), 100 ml capacity Wick: Includes one, 0.4 cm x 15 cm (0.15625 in x 5.875 in)	581148

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE WUFI OR PRODUCT INSERT.

Questions: bvimedical.com/customer-support



BVI
bvimedical.com

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2023

1721526-01 IN