

a new era has just begun



ONE-STEPYOGURT TECHNIQUE

EASY

TREPHINATION EXTENDS UP TO SCHWALBE'S LINE IN THE PERIPHERY WHERE DM IS VERY EASILY PEELED OFF.

THE TECHNIQUE RESEMBLES, THUS, THE
OPENING OF A **YOGURT** CUP CREATING A HINGE
TO GRASP THE GRAFT.

SAFE

NO ENDOTHELIAL CELL LOSS AS ALL MANIPULATIONS ARE MADE ONLY ON EXCESS TISSUE WHICH ALSO ACTS AS AN ORIENTATION MARKING.

NO MORE UPSIDE-DOWN
PLACED GRAFT

PATENT PENDING

ONE-STEP
YOGURT TECHNIQUE

FIXED DEPTH

available size

JD 700 TZ Ø 7.0

JD 750 TZ Ø 7.5

JD 800 TZ Ø 8.0

JD 850 TZ Ø 8.5

JD 900 TZ Ø 9.0

OUICK

PREPARATION IS PERFORMED IN A
SINGLE-PEEL TECHNIQUE IN LESS THAN
5 MINUTES

AFFORDABLE

NO NEED FOR SPECIAL INSTRUMENTATION OR MULTIPLE CORNEAL PUNCHES

DO IT YOUR SELF

CREATE YOUR GRAFT IN OPERATING THEATRE
VERY LOW LEARNING CURVE







Sistema di Qualità Certificato Certified Quality System

UNI EN ISO 9001 UNI CEI EN ISO 13485 © Copyright: Hermann Janach



Via Borgo Vico, 35 - 22100 Como - Italy Tel. +39 031 574088 - fax +39 031 572055 international@janach.it

BRO.086.07.19.IT/EN