

DualBore SideFlō™ Cannula

NEW



- New design **increases safety** by using a side port for fluid injection, eliminating possible retinal injury caused by a fluid jet stream.
- Multiple fluid egress vents allow double the outflow over the previous model for **faster pressure relief** during injection.
- All-metal exterior allows **easy entry through valved cannula** systems.
- Special hub decreases air bubbles during PFCL injection so you can **complete the case without interruption**.
- Optimized for injection of PFCL, staining dyes, or other liquids in small-gauge vitrectomy surgery.

3423 DualBore SideFlō™ Cannula 23g
3425 DualBore SideFlō™ Cannula 25g
3427 DualBore SideFlō™ Cannula 27g

Patent Pending

Contact us today for more information or to arrange an evaluation!



See it in action!
Scan the code or visit us at
MedOne.com to see a video

MedOne Surgical, Inc.
670 Tallevast Road
Sarasota, FL 34243

941.359.3129 (TEL)
866.633.6631 (USA)
www.MedOne.com
MedOne@MedOne.com

MedOne
Exactly what you're looking for