

Backflush Handle

Our Backflush Handle is designed to maximize performance and ergonomic comfort. The short handle provides excellent maneuverability in the surgical field, especially around wide-field viewing systems. The grip area is designed for precise control in delicate procedures. The luer lock fitting enables a quick, easy connection with any luer-compatible cannula, including all MedOne cannulas. The removable plug allows use with active or passive aspiration.

3228 Backflush Handle



Backflush Instruments

Backflush Instruments combine our ergonomic handle with our most popular cannulas in one convenient instrument. The instruments are available with a blunt tip cannula or a .75mm silicone tip, allowing easy passage through valved entry cannulas. The removable plug allows use with active or passive aspiration.

3286 Backflush Cannula 23g

3290 Backflush FlexTip™ 23g

3295 Backflush Cannula 25g

3296 Backflush FlexTip™ 25g

3276 Backflush Cannula 27g

3278 Backflush FlexTip™ 27g

