













Beaver Mini-Blades

Available in a variety of designs to meet all small incision and excision needs. See page 26-27 for Beaver Handles and Accessories.

| | |
|--|---|
| 374562 Mini-Blade Orthopedic Mini-Blade, Length 1.255 inches (31.88 mm).. 24 per box. |  |
| 376100 Mini-Blade Overall length 1.31 inches (33.27 mm). 24 per box. |  |
| 376300 Mini-Blade Sharp on two sides, with pointed tip. 6 per box. |  |
| 376401 Mini-Blade MIS 6400. Narrower version of Beaver6400. For minimally invasive surgery. 24 per box. |  |
| 376500 Mini-Blade Pointed tip, smaller alternative to #11 general blade. 24 per box. |  |
| 374769 Mini-Blade Knife Surgical blade, 180 degree cutting surface. 3 per box. |  |

Grindless Beaver Mini-Blades*

Beaver-Visitec's Mini-Blades are manufactured using a proprietary grindless technology that yields consistent sharpness and blade geometry and are available in a variety of designs to meet all small incision and excision needs. Beaver Mini-Blades are preassembled into a unique, patented tip, designed for safe assembly and ease of use. Each Mini-Blade is color coded for easy identification. See page 26-27 for Beaver Handles and Accessories.

| | |
|---|--|
| Beaver6400 6400 Mini-Blade Grindless, rounded tip, sharp on one side, green. 20 per box. |  |
| Beaver6600 6600 Mini-Blade Grindless, lamellar, angled 60 degrees, smooth undersurface, sharp all around, bevel up, black. 3 per box. |  |
| Beaver6610 6610 Mini-Blade Grindless, lamellar, straight, smooth undersurface, sharp all around, black. 20 per box. |  |
| Beaver6700 6700 Mini-Blade Grindless, smaller alternative to #15 general blade, orange. 20 per box. |  |
| Beaver6900 6900 Mini-Blade Grindless, straight, sharp all around, 180 degree cutting surface, double bevel blade, blue. 20 per box. |  |
| Beaver6910 6910 Mini-Blade Grindless, cleft palate blade, angled 60 degrees, sharp all around, 180 degree cutting surface, double bevel blade, blue. 6 per box. |  |

* Patent D636082. Additional US and International patents pending.