



# Needle Holders



# Needle Holders

72051

## Castroviejo Needle Holder




8 mm straight jaws. Without lock.  
Overall length 142 mm.



72052

## Castroviejo Needle Holder



8 mm straight jaws. With lock.  
Overall length 142 mm.



72053

## Castroviejo Needle Holder



8 mm curved jaws. Without lock.  
Overall length 142 mm.



72054

## Castroviejo Needle Holder



8 mm curved jaws. With lock.  
Overall length 142 mm.



73021

## Barraquer Needle Holder



Fine 9 mm Gently curved jaws. Without lock.  
Overall length 134 mm.



73022

## Barraquer Needle Holder



Fine 9 mm Gently curved jaws. With lock.  
Overall length 134 mm.



73023

## Barraquer Needle Holder



Fine 9 mm straight jaws. Without lock.  
Overall length 136 mm.



73030

## Barraquer Needle Holder



Fine 9 mm Gently curved jaws. Without lock. Overall length 135 mm.



73050

## Barraquer Needle Holder



9 mm Gently curved jaws. Without lock. Overall length 116 mm.



73111

## Barraquer Needle Holder



Fine 10 mm Gently curved jaws. Inner surfaces roughened. Handle with lengthwise slots and screwed end-joint. Without lock. Overall length 125 mm.



3-205

## STORZ Needle Holder 7 mm Curved Fine



Curved 7 mm fine jaws. Round knurled handle with dull finish. Titanium. Overall length 106 mm.



33-203

## Needle Holder 9 mm Curved Medium Long Handle



Curved 9 mm medium jaws. Round knurled handle with dull finish. Without lock. Titanium. Overall length 137 mm.



E3796

## STORZ Needle Holder 9 mm Straight Medium



Medium 9 mm jaws. Large round knurled handle. Polished finish. Without lock. Overall length 130 mm.



E3807 WO

## Osher Needle Holder



Delicate curved and tapered 9 mm jaws for grasping fine 10-0 sutures. Without lock. Overall length 110 mm.



# Needle Holders

E3828 CW

## Barraquer Needle Holder



Delicate 9 mm Gently curved jaws. With lock.  
Overall length 133 mm.



E3828 CWO

## Barraquer Needle Holder



Delicate 9 mm Gently curved jaws. Without lock.  
Overall length 133 mm.



E3828 SW

## Barraquer Needle Holder



Delicate 9 mm straight jaws. With lock.  
Overall length 134 mm.



E3828 SWO

## Barraquer Needle Holder



Delicate 9 mm straight jaws. Without lock.  
Overall length 134 mm.



E3835

## McPherson Needle Holder



Fine 10 mm moderately curved jaws. Large handle.  
Fine 11 mm strongly curved jaws. With lock.  
Overall length 107 mm.



E3836 C

## Cohan Needle Holder



Fine 8 mm Gently curved jaws. Without lock.  
Overall length 104 mm.



E3836 CW

## Cohan Needle Holder



Fine 8 mm Gently curved jaws. With lock.  
Overall length 104 mm.



E3839

## Troutman Needle Holder 10 mm Curved Medium



Medium 10 mm gently curved jaws. Medium round horizontally serrated handle with dull finish. Without lock. Overall length 122 mm.



E3842

## Barraquer Needle Holder



Fine 10 mm moderately curved jaws. Large handle. Without lock. Overall length 137 mm.



E3843

## Barraquer Needle Holder



Fine 10 mm moderately curved jaws. Large handle. With lock. Overall length 137 mm.



E3844

## Barraquer Needle Holder



Fine 10 mm Gently curved jaws. Medium size handle. Without lock. Overall length 135 mm.



E3845

## Barraquer Needle Holder



Fine 10 mm Gently curved jaws. wide short handle. Without lock. Overall length 125 mm.



E3848

## Castroviejo Needle Holder



Medium 11 mm Gently curved jaws. With lock. Overall length 137 mm.



E3849

## Castroviejo Needle Holder



Medium 11 mm Gently curved jaws. Without lock. Overall length 137 mm.



# Needle Holders

E3850

## Castroviejo Needle Holder



Heavy 12 mm straight jaws. With lock.  
Overall length 138 mm.



E3852

## Castroviejo Needle Holder



Medium 9 mm tungsten carbide inlaid straight jaws.  
With lock. Overall length 148 mm.



E3861

## Castroviejo Needle Holder



Medium 10 mm straight jaws. With lock.  
Overall length 136 mm.



E3861 WHMC

## Castroviejo Needle Holder



Fine 9 mm curved jaws. With lock.  
Overall length 130 mm.



E3864

## Castroviejo Needle Holder



Heavy 16 mm straight jaws. With lock.  
Overall length 149 mm.



E3876

## Derf Needle Holder



Heavy 13 mm serrated straight jaws with grooves.  
Ring handle with ratchet lock. Overall length 130 mm.

