


# Oil Removal Simplified

Quick & easy oil removal with specialized trocar cannula

**NEW**



## 3300 Oil Removal Cannula (23g) Kapran

Simply place the trocar cannula using standard technique, then push the included silicone tube onto the specially shaped port cannula. Attach a syringe and aspirate the oil using a standard syringe or pump system. Enables removal of oil through a sutureless incision.

Developed in cooperation with Ziya Kapran, MD, Istanbul, Turkey

**Contact us today** for more information or to arrange an evaluation!

**MedOne Surgical, Inc.**  
670 Tallevast Road  
Sarasota, FL 34243

941.359.3129 (TEL)  
866.633.6631 (USA)  
[www.MedOne.com](http://www.MedOne.com)

[MedOne@MedOne.com](mailto:MedOne@MedOne.com)



**See it in action!**  
Scan the code or visit us at  
[MedOne.com](http://MedOne.com) to see a video

 **MedOne**  
*Exactly what you're looking for*