

Green PolyTip[®] Cannula

NEW


Green polyimide tip allows enhanced visibility for accurate positioning during surgery

Short 2mm tip length results in increased stiffness for easier insertion into the subretinal space

Available in 23g, 25g, and 27g



- 3254** PolyTip[®] Cannula 23g/38g (2mm)
- 3255** PolyTip[®] Cannula 25g/38g (2mm)
- 3256** PolyTip[®] Cannula 27g/38g (2mm)



MedOne Surgical, Inc.
670 Tallevast Road
Sarasota, FL 34243

941.359.3129 (TEL)
866.633.6631 (USA)
www.MedOne.com
MedOne@MedOne.com

*Call now to order or
to arrange an evaluation!*

 **MedOne**[®]
Exactly what you're looking for