



Retractors

- Iris 156
- Lacrimal..... 156
- Lid..... 157
- Orbital..... 158



Retractors

Iris

54055

Iris Retractor With Flexible Fixing

Silicone slide for flexible fixation of the iris. 7 mm.



x 1

E0750 NS

MADLAB Iris Retractor 5/pack

Developed for ease of use, these very flexible retractors are made of the super-elastic Nitinol alloy to provide durability and consistency. The retractors are supplied in packs of five, are fully reusable and are shipped in a unique sterilization case. Overall length 10 mm.



x 5

Lacrimal

57061

Axenfeld Lacrimal Sac Retractor

BAUSCH + LOMB
Instruments

Blunt. Overall length 142 mm.



x 1

57071

Rollet Lacrimal Sac Retractor

BAUSCH + LOMB
Instruments

Blunt. Overall length 120 mm.



x 1

E4500

Stevenson Lacrimal Sac Retractor

BAUSCH + LOMB
Instruments

3 mm teeth. Blade spread 20 mm.
Overall length 85 mm.



x 1

Lid

52051

Desmarres Lid Retractor

BAUSCH + LOMB
Instruments



Size 1. 11 mm blade. Overall length 140 mm.



52191

Desmarres Lid Retractor

BAUSCH + LOMB
Instruments



Baby Size. 9 mm blade. Overall length 140 mm.



E0980

Desmarres Lid Retractor

BAUSCH + LOMB
Instruments



Size 0. 11 mm. Overall length 140 mm.



E0981

Desmarres Lid Retractor

BAUSCH + LOMB
Instruments



Size 1. 13 mm. Overall length 137 mm.



E0982

Desmarres Lid Retractor

BAUSCH + LOMB
Instruments



Size 2. 15 mm. Overall length 143 mm.



E0997

Jaffe Lid Retractor

BAUSCH + LOMB
Instruments



15 mm blades, one pair. Overall length 53 mm.



E1064

Jaeger Lid Plate

BAUSCH + LOMB
Instruments



20 mm x 35 mm and 24 mm x 35 mm curved ends.
Clear lucite. Overall length 93 mm.



Retractors

Lid

E1069

Berke-Jaeger Lid Plate

BAUSCH + LOMB
Instruments



Curved over length and width.
20 mm and 23 mm wide blades. Metal with dull finish.
Overall length 109 mm.



E1069S

Khan-Jaeger Lid Plate

BAUSCH + LOMB
Instruments



Designed for transconjunctival blepharoplasty.
Non-reflective surface. 20 mm and 23 mm wide blades.
Overall length 109 mm.



Orbital

E3705

Orbital Retractor-Elevator

BAUSCH + LOMB
Instruments



Flexible. Double-ended 16 + 25 mm ends.
Overall length 182 mm.



E3706

Orbital Retractor-Elevator

BAUSCH + LOMB
Instruments



Flexible. Double-ended 10 + 14 mm ends.
Overall length 180 mm.



E5008

Schepens Orbital Retractor

BAUSCH + LOMB
Instruments



14 x 56 mm blade. Overall length 148 mm.

