

Freedom of choice to optimize your surgery

ARCALINE

98% Perfluorodecaline

Due to its specifications, Arcaline is particularly recommended in association with tamponade of silicone oil

ARCOANE

99% Perfluoro-n-octane

Due to its specifications, Arcotane is particularly recommended in association with tamponade of gases

Comments	Advantages	Specifications		Specifications	Advantages	Comments
Stops bleeding and stabilizes the retina	Stronger pressure	1,94	Density	1,79	Efficient	In simple detachment
Photo coagulation aid	Few bubbles	141°C	Boiling point	110°C		
		1,31 at 20°C	Refractive index	1,29 at 20°C	Better visualization	Easy to identify
Post surgery residual bubble	Leaves a film on the retina	2,4 cst	Viscosity	0,8 cst	Rapidly drops on to the bottom of the eye	Allows complete withdrawal
Post surgery residual bubble	0,1 g evaporates within 15 min	13 mm Hg	Vapor pressure	56 mm Hg	0.1g evaporates within 2 min	Disappears during Pfc/air exchange

Manufactured in accordance with Pharmaceutical industry standards for injectable solutions.

	FDA standard PfcI	ARCAD STANDARD
Fluoride	≤ 2ppm	No detectable
Non volatile Residue	≤ 20ppm	No detectable
C-H link	≤ 0,015	No detectable

STERILE - SINGLE USE - APYROGENIC

ARCALINE



ARCOANE



ARCAD
O P H T A®

Arcadophta 11, rue Antoine Ricord
31100 TOULOUSE - FRANCE
Tél. : +33 5 61 40 52 35
Fax : +33 5 61 40 84 66
E-mail : info@arcadophta.com
Skype : arcad4

Find updated bibliographies on the website www.arcadophta.com