


## Trays

### Polymer – removable lid

#### 0112013

Tray holds 10 instruments with silicone finger matting. Outer size W260mm D160mm H25mm. Inner size W250mm D150mm H20mm.




#### 0112015

Tray holds 20 instruments with silicone finger matting. Outer size W390mm D265mm H25mm. Inner size W380mm D250mm H20mm.




#### 0112014

Tray, 2 tier, holds 20 instruments with silicone finger matting. Outer size W260mm D160mm H45mm. Inner size W250mm D150mm H40mm.



#### 0112019

Tray, 2 tier, holds 40 instruments with silicone finger matting. Outer size W390mm D270mm H45mm. Inner size W380mm D255mm H40mm.




# Trays

Polymer – removable lid



**0112012**

Diamond Knife tray with silicone grip for 1 knife. Outer size W160mm D70mm H25mm. Inner size W150mm D60mm H20mm.



**0112016**

Tray holds 1 - 2 instruments with silicone finger matting. Outer size W160mm D70mm H25mm. Inner size W150mm D60mm H20mm.



**0112017**

Tray holds 3 instruments with silicone finger matting. Outer size W175mm D110mm H25mm. Inner size W160mm D100mm H20mm.



**0112018**

Tray holds 3 instruments with silicone finger matting. Outer size W200mm D110mm H25mm. Inner size W190mm D100mm H20mm.

Number of instruments per tray is a guide only.

## Trays

Metal – removable lid

### 0112003

Tray holds 12 instruments with 3 strips of silicone finger matting. Outer size W220mm D210mm H30mm. Inner size W200mm D200mm H25mm.



### 0112002

Tray holds 20 instruments with 4 strips of silicone finger matting. Outer size W385mm D270mm H30mm Inner size W365mm D260mm H25mm.



### 0112004

Tray with 2 areas. Area 1 holds 20 instruments with 4 strips of silicone finger matting.

Size W280mm D260mm H25mm.

Area 2 With silicone finger matting accommodates larger items such as a phaco handpiece, size W260mm D160mm H50mm.

Total size of tray outer W510mm D270mm H70mm. Inner W500mm D260mm H65mm.



Number of instruments per tray is a guide only.