



Trephines/ Punches

LASEK	208
Miscellaneous	208
Punches	209
Rongeurs	209
Trephines	209



Trephines/Punches

LASEK

53008T

Trephine LASEK Titanium, 8.0 mm



The trephine has a 270° cutting edge. A cut-out of 90° allows to create a flap hinge. Marks on top of the ring indicate the hinge area. Titanium 8.0 mm Ø.



53009T

Trephine LASEK Titanium, 9.0 mm



The trephine has a 270° cutting edge. A cut-out of 90° allows to create a flap hinge. Marks on top of the ring indicate the hinge area. Titanium 9.0 mm Ø.



Miscellaneous

E3048

Fine Corneal Carrying Case



Plastic case equipped with slots for arranging sutures.



E3067

Brightbill Corneal Cutting Block



Incorporates three wells of various sizes with red 8 mm centering dot: Small, 15 mm diameter; Medium, 16.5 mm diameter; Large, 18.0 mm diameter. Wells on each side for added instrument life. Size 52 x 52 x 19 mm.



93100S

Böhnke Donor Cornea Holder



The Böhnke Donor Cornea Holder is used for holding the donor cornea during cornea storage in cell culture medium. Packaged single. Sterile. 10/box.



Punches

E2798

Kelly Descemet's Membrane Punch



For 0.75 mm round holes without tissue tags. Triangular squeeze action handle. Overall length 130 mm.



E3097 BP

Base Plug



For Universal Trephine Handle with Cottingham Punch (E3097). 14 mm in diameter, 8 mm height.



E3097

Universal Trephine Handle with Cottingham Punch



Metal punch with universal style handle. Uses E3096 Series STORZ Corneal Transplant Trephine Blades. Supplied with two replaceable base plugs (E3097 BP). Dull finish. Overall length 60 mm.



Rongeurs

E4590

Belz Lacrimal Sac Rongeur

BAUSCH + LOMB
Instruments



Rounded 17 mm long jaws. 2 mm wide at tip, 3 mm wide at midline. Overall length 180 mm.



Trephines

E3050 6.0G

Tanne Universal Disposable Trephine



6.0 mm Ø. Sterile. For handle E3050 H.



Trephines/Punches

Trephines

E3050 6.25G

Tanne Universal Disposable Trephine

6.25 mm Ø. Sterile. For handle E3050 H.



E3050 6.5G

Tanne Universal Disposable Trephine

6.5 mm Ø. Sterile. For handle E3050 H.



E3050 6.75G

Tanne Universal Disposable Trephine

6.75 mm Ø. Sterile. For handle E3050 H.



E3050 7.0G

Tanne Universal Disposable Trephine

7.0 mm Ø. Sterile. For handle E3050 H.



E3050 7.25G

Tanne Universal Disposable Trephine

7.25 mm Ø. Sterile. For handle E3050 H.



E3050 7.5G

Tanne Universal Disposable Trephine

7.5 mm Ø. Sterile. For handle E3050 H.



E3050 7.75G

Tanne Universal Disposable Trephine

7.75 mm Ø. Sterile. For handle E3050 H.



Trephines

E3050 8.0G

Tanne Universal Disposable Trephine

8.0 mm Ø. Sterile. For handle E3050 H.



E3050 8.25G

Tanne Universal Disposable Trephine

8.25 mm Ø. Sterile. For handle E3050 H.



E3050 8.5G

Tanne Universal Disposable Trephine

8.5 mm Ø. Sterile. For handle E3050 H.



E3050 8.75G

Tanne Universal Disposable Trephine

8.75 mm Ø. Sterile. For handle E3050 H.



E3050 9.0G

Tanne Universal Disposable Trephine

9.0 mm Ø. Sterile. For handle E3050 H.



E3050 9.25G

Tanne Universal Disposable Trephine

9.25 mm Ø. Sterile. For handle E3050 H.



E3050 H

Universal Handle for Blades

Designed for use with the E3050 series Tanne Disposable Trephine Blades. Overall length 60 mm.



Trephines/Punches

Trephines

E3095

STORZ Universal Trephine Handle

Round knurled handle with adjustable finger-like tips that rest in the grooved disposable trephine blades. Used with E3096 Series Disposable Trephine Blades. Polished finish. Overall length 49 mm.



E3096 6.0G

Disposable Trephine Blade

Sterile. 6.0 mm Ø. Length 9.0 mm.
To be used with handle E3095.



E3096 6.0LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
6.0 mm Ø. 9.0 mm.



E3096 6.2G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
6.25 mm Ø. 9.0 mm.



E3096 6.5G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
6.5 mm Ø. 9.0 mm.



E3096 6.5LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
6.5 mm Ø. 15.0 mm.



E3096 6.7G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
6.75 mm Ø. 9.0 mm.



Trephines

E3096 6.7LG	<p>Disposable Trephine Blade</p> <p>Long. Sterile. For holder E3095 or punch E3097. 6.75 mm Ø, 15.0 mm.</p>		 
E3096 7.0G	<p>Disposable Trephine Blade</p> <p>Sterile. For holder E3095 or punch E3097. 7.0 mm Ø, 9.0 mm.</p>		 
E3096 7.0LG	<p>Disposable Trephine Blade</p> <p>Long. Sterile. For holder E3095 or punch E3097. 7.0 mm Ø, 15.0 mm.</p>		 
E3096 7.2G	<p>Disposable Trephine Blade</p> <p>Sterile. For holder E3095 or punch E3097. 7.25 mm Ø, 9.0 mm.</p>		 
E3096 7.2LG	<p>Disposable Trephine Blade</p> <p>Long. Sterile. For holder E3095 or punch E3097. 7.25 mm Ø, 15.0 mm.</p>		 
E3096 7.5G	<p>Disposable Trephine Blade</p> <p>Sterile. For holder E3095 or punch E3097. 7.5 mm Ø, 9.0 mm.</p>		 
E3096 7.5LG	<p>Disposable Trephine Blade</p> <p>Long. Sterile. For holder E3095 or punch E3097. 7.5 mm Ø, 15.0 mm.</p>		 

Trephines/Punches

Trephines

E3096 7.7G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
7.75 mm Ø. 9.0 mm.



E3096 7.7LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
7.75 mm Ø. 15.0 mm.



E3096 8.0G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
8.0 mm Ø. 9.0 mm.



E3096 8.0LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
8.0 mm Ø. 15.0 mm.



E3096 8.2G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
8.25 mm Ø. 9.0 mm.



E3096 8.2LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
8.25 mm Ø. 15.0 mm.



E3096 8.5G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
8.5 mm Ø. 9.0 mm.



Trephines

E3096 8.5LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
8.5 mm Ø. 15.0 mm.



E3096 8.7G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
8.75 mm Ø. 9.0 mm.



E3096 8.7LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
8.75 mm Ø. 15.0 mm.



E3096 9.0G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
9.0 mm Ø. 9.0 mm.



E3096 9.0LG

Disposable Trephine Blade

Long. Sterile. For holder E3095 or punch E3097.
9.0 mm Ø. 15.0 mm.



E3096 9.2G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
9.25 mm Ø. 9.0 mm.



E3096 9.5G

Disposable Trephine Blade

Sterile. For holder E3095 or punch E3097.
9.5 mm Ø. 9.0 mm.



