

auroSphere



For Enucleation & Evisceration procedure

- Restores the natural appearance of the eye and surrounding tissues
- Enables both motility and cosmesis with better accommodation
- Silicone material of the sphere makes it easy to suture with extra – ocular muscles.
- Aurosphere box packed with one sphere ball and two foils of 6-0 polycryl (absorbable) suture

Specification:

Product Code	AS - 14	AS - 16	AS - 18	AS - 20	AS - 22
Size (Dia)	14 mm	16 mm	18 mm	20 mm	22 mm