



Scissors

Capsule	160
Conjunctival	161
Corneal	162
Enucleation	167
Foreign Body	167
Intraocular Lens	168
Iris	168
Miscellaneous	169
Retinal	172
Strabismus	173
Suture Stich	174
Tenotomy	175



Scissors


Capsule

24081

Vannas Capsulotomy Scissors, Straight



Overall length 82 mm.




24082

Vannas Capsulotomy Scissors, Curved



85 mm.




24083

Vannas Capsulotomy Scissors, Straight



Model Tuebingen. Overall length 86 mm.



24084

Vannas Capsulotomy Scissors, Curved



Model Tuebingen. Overall length 85 mm.



25030

Ong Capsule Scissors



Angled up. 11 mm from mid-screw to tip. 85 mm.



25040

Uthoff Capsule Scissors



Curved. 83 mm.



25030Y2

Ong Capsule Scissors




Slightly curved tip. Angled.



Capsule

E3386

Vannas Straight Capsulotomy Scissors




Sharp pointed tips. 8 mm from mid-screw to tip.
Overall length 84 mm.



E3387

Vannas Curved Capsulotomy Scissors



6mm blades with sharp pointed tips. 8 mm from mid-screw to tip. Serrated handle with polished finish.
Overall length 84 mm.



E3388

Vannas Angled Capsulotomy Scissors



6mm blades with sharp pointed tips. 9 mm from mid-screw to tip. Serrated handle with polished finish.
Overall length 81 mm.



E3388 GWL

Gills-Welsh Vannas Micro Scissors



Angled. Long model. 11 mm blades with sharp pointed tip. 13 mm from mid-screw to tip. Serrated handle with polished finish. Overall length 86 mm.



E3390 A

Angled Capsule Scissors



Delicate 10 mm blades angled 45°. 110 mm.



Conjunctival

E3315

Troutman Curved Conjunctival Scissors



Blunt tips. 17 mm from mid-screw to tip. With lock.
Overall length 112 mm.



Scissors

Conjunctival

E3316

Westcott Conjunctival Scissors

Lightly curved. Blunt tips. 16 mm from mid-screw to tip.
Overall length 132 mm.



E3450

Wilmer Conjunctival and Utility Scissors

Angled sharp pointed tips. 18 mm from mid-screw to tip.
Overall length 94 mm.



Corneal

23010

Castroviejo Curved Corneal Scissors


Blunt tips. 18 mm from mid-screw to tip.
Overall length 94 mm.



23121

Troutman Microscissors


Right model. Blunt tips. Overall length 99 mm.



23122

Troutman Microscissors


Left model. Blunt tips. Overall length 99 mm.



23130

Micro Scissors

Blade length 10 mm. Cutting length 8 mm. pointed tips.
Straight. Overall length 100 mm.



Corneal

23140

Micro Scissors




Sharp/sharp. Slightly curved. Length of cut 6 mm. Overall length 105 mm.



23160

Micro Scissors



Angled up. Overall length 100 mm.



23201

Katzin Corneal Transplant Scissors



Right. Overall length 98 mm.



23202

Katzin Corneal Transplant Scissors



Left. Overall length 98 mm.



E3217 L

Fogla Scissors, Left



These scissors are used to remove the residual corneal stroma separated from the Descemet's membrane in DALK surgery. The lower blade of the scissors is thinner and 1 mm longer than the upper blade which displaces Descemet's membrane, protecting it from being damaged by the blades. Left. Overall length 107 mm.



E3217 R

Fogla Scissors, Right



This scissor is used to remove the residual corneal stroma separated from the Descemet's membrane in DALK surgery. The lower blade of the scissors is thinner and 1 mm longer than the upper blade which displaces Descemet's membrane, protecting it from being damaged by the blades. Right. Overall length 107 mm.



Scissors

Corneal

E3218 L

Castroviejo Corneal Scissors, Left




Mini. Curved. Blunt tips. 11 mm from mid-screw to tip.
Overall length 106 mm.



E3218 MM

Castroviejo Mini Corneal Scissors



Universal. Blunt tips. 11 mm from mid-screw to tip.
Overall length 106 mm.



E3218 R

Castroviejo Corneal Scissors, Right



Mini. Curved. Blunt tips. 11 mm from mid-screw to tip.
Overall length 106 mm.



E3220

Castroviejo Curved Corneal Scissors



Blunt tips. 16 mm from mid-screw to tip.
Overall length 115 mm.



E3221

Castroviejo Keratoplasty Scissors



Angled. Blunt tips. Lower blade 1 mm longer. 16 mm from
mid-screw to tip. Overall length 110 mm.



E3222 R

Castroviejo Corneal Scissors



Right orientation, knee angle. Blunt tips. 16 mm from
mid-screw to tip. Overall length 106 mm.



E3223

Castroviejo Angled Corneal Scissors



Universal mini model. Blunt tips. 19 mm from mid-screw
to tip. Overall length 106 mm.



Corneal

E3231L

Mini Corneal Scissors, Left



Delicate. For completion of trepanations.
Lower blade 0.5 mm longer. 9 mm from mid-screw to tip.
Overall length 100 mm.



E3231R

Mini Corneal Scissors, Right



Delicate. For completion of trepanations.
Lower blade 0.5 mm longer. 9 mm from mid-screw to tip.
Overall length 100 mm.



E3232

Katzin Corneal Transplant Scissors



Left, miniature model. Strongly curved. Blunt tips. Lower blade 12 mm, upper blade 11 mm (from mid-screw to tip).
Overall length 108 mm.



E3233

Katzin Corneal Transplant Scissors



Right, miniature model. Strongly curved. Blunt tips. Lower blade 12 mm, upper blade 11 mm (from mid-screw to tip).
Overall length 108 mm.



E3300

Castroviejo Corneal Section Scissors



Right. Angled. Blunt tips. Curved blades. Lower blade 1 mm longer and rounded spatula-like to minimise trauma. 18 mm from mid-screw to tip.
Overall length 110 mm.



E3301

Castroviejo Corneal Section Scissors



Left. Angled. Blunt tips. Curved blades. Lower blade 1 mm longer and rounded spatula-like to minimise trauma. 18 mm from mid-screw to tip.
Overall length 110 mm.



E3304

Castroviejo Corneal Section Scissors



Right. Angled. Blunt tips. Curved blades. Lower blade 0.5 mm long, 14 mm from mid-screw to tip.
Overall length 106 mm.



Scissors

Corneal

E3304M

Castroviejo Corneal Section Scissors



Right. Curved. Blunt tips. Lower blade 0.5 mm longer. 13 mm from mid-screw to tip. Overall length 106 mm.



E3305

Castroviejo Corneal Section Scissors



Left. Angled. Blunt tips. Curved blades. Lower blade 0.5 mm long. 14 mm from mid-screw to tip. Overall length 106 mm.



E3305M

Castroviejo Corneal Section Scissors



Left, mid-mini. Curved. Blunt tips. Lower blade 0.5 mm longer. 13 mm from mid-screw to tip. 106 mm.



E3306

Troutman-Castroviejo Corneal Scissors



Mini model. Curved. Blunt tips. 13 mm from mid-screw to tip. Handle with lock. Overall length 116 mm.



E3307

Troutman-Castroviejo Corneal Scissors



Mini model. Curved. Blunt tips. 13 mm from mid-screw to tip. Handle with lock. Overall length 116 mm.



E3321MWH

Castroviejo Keratoplasty-Utility-Stitch Scissors



Very small pointed tips. 13 mm from mid-screw to tip. Overall length 131 mm.



E3490

Aebli Corneal Scissors



Straight. Blunt tips. Overall length 92 mm.



Enucleation

E3651

Curved Enucleation Scissors

Medium rounded blades. Blunt tips. 44 mm from mid-screw to tip. Overall length 135 mm.



E3652

Curved Enucleation Scissors

Heavy rounded blades. Blunt tips. 44 mm from mid-screw to tip. Overall length 130 mm.



Foreign Body

ET8306CTL

Snyder AC Scissors Tip Curved



These scissors offer control through a small incision and can be used for cutting fine structures including capsule, iris, adhesions and sutures in the anterior chamber. Fine tapered and gently curved blades provide precise cutting and can be rotated left or right for enhanced control in the eye. The 23 Gauge shaft diameter is ideal for entry through a small stab incision or a phaco wound.

Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 49 mm.



ET8306STL

Snyder AC Scissors Tip Straight



These scissors offer control through a small incision and can be used for cutting fine structures including capsule, iris, adhesions and sutures in the anterior chamber. Fine tapered blades provide precise cutting while maximizing visibility in the eye. The 23 Gauge shaft diameter is ideal for entry through a small stab incision or a phaco wound. Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 49 mm.



Intraocular Lens

ET1306

Osher IOL Scissors angled



The blades of this scissors are narrow to provide ease of entry through the incision. Serrations in the blade allow secure cutting of the implant without excessive movement. Ideal for acrylic and silicone implants. Made from titanium. Overall length 91 mm.



Iris

24071

Noyes Iris Scissors



Blunt/blunt. Overall length 113 mm.



E3343

Straight Iris Scissors



Sharp pointed tips, 21 mm from mid-screw to tip. Overall length 121 mm.



E3344

Curved Iris Scissors



Sharp pointed tips, 21 mm from mid-screw to tip. Overall length 121 mm.



E3347

McPherson-Vannas Curved Iris Scissors



Sharp pointed tips, 11 mm from mid-screw to tip. Overall length 83 mm.



E3349

Cohan-Vannas Curved Iris Scissors



Sharp pointed tips, 12 mm from mid-screw to tip. Overall length 57 mm.



Iris

E3360

DeWecker Iris Scissors

Blunt/sharp pointed tips. 12 mm from bend to tip.
Overall length 114 mm.



E3366

Barraquer Iris Scissors

Blunt/blunt. 8 mm from bend to tip.
Overall length 54 mm.



E3404

Straight Iris Scissors

Sharp pointed tips. 28 mm from mid-screw to tip.
Overall length 100 mm.



E3405

Curved Iris Scissors

Sharp pointed tips. 28 mm from mid-screw to tip.
Overall length 100 mm.



E3408

Curved Iris Scissors

Light standard pattern. Sharp pointed tips. 25 mm from mid-screw to tip. Overall length 98 mm.



Miscellaneous

23060

Bonn Eye Scissors

Curved. Overall length 89 mm.



Scissors

Miscellaneous

24020

Eye Scissors, curved

Delicate. Overall length 105 mm.



24050

Eye Scissors

Curved. Overall length 90 mm.



E3243

Clayman-Westcott Curved Scissors

Blunt tips. 11 mm from mid-screw to tip.
Overall length 107 mm.



E3322

Westcott Type Utility Scissors

Rounded sharp tips. 25 mm from mid-screw to tip.
Overall length 131 mm.



E3383 C

Clayman-Vannas Curved Scissors

Sharp pointed tips. 7 mm from mid-screw to tip.
Overall length 83 mm.



E3383 S

Clayman-Vannas Straight Scissors

Sharp pointed tips. 7 mm from mid-screw to tip.
Overall length 83 mm.



Miscellaneous

E34101

Straight Eye Scissors

Pointed tips. 29 mm from mid-screw to tip.
Overall length 109 mm.



E34121

Curved Eye Scissors

Pointed tips. 29 mm from mid-screw to tip.
Overall length 109 mm.



E3424

Straight Eye Scissors

Pointed tips. 32 mm from mid-screw to tip.
Overall length 109 mm.



E3426

Curved Eye Scissors

Pointed tips. 32 mm from mid-screw to tip.
Overall length 109 mm.



E3436

Angled Eye Scissors

Blunt probe-point tips. Angled. 26 mm from mid-screw to tip.
Overall length 116 mm.



E3580

Littler Curved Dissecting Scissors

Light curved. Blunt tips. eyes in blades used for drawing sutures or muscle through a tunnel dissection.
27 mm from mid-screw to tip. Overall length 119 mm.



Retinal

ET8315

23 Gauge StableGrip™ Vertical Scissors



The sharp pointed rotatable tips of this scissors are designed to facilitate removal of membranes that are causing traction. These scissors are designed for PVR and tangential traction. The front half of this instrument is gold to indicate a vertical scissors. The back half is colored green to indicate 23 Gauge. Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 58 mm.



ET8316

23 Gauge StableGrip™ Curved Scissors



These scissors are designed for delamination of epiretinal tissue to release traction. Can also be used for PVR. The front half of this instrument is pink to indicate a curved scissors. The back half is colored green to indicate 23 Gauge. Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 56 mm.



ET8515

25 Gauge StableGrip™ Vertical Scissors



The sharp pointed rotatable tips of this scissors are designed to facilitate removal of membranes that are causing traction. These scissors are designed for PVR and tangential traction. The front half of this instrument is gold to indicate a vertical scissors. The back half is colored dark blue to indicate 25 Gauge. Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 58 mm.



ET8516

25 Gauge StableGrip™ Curved Scissors



These scissors are designed for delamination of epiretinal tissue to release traction. Can also be used for PVR. The front half of this instrument is pink to indicate a curved scissors. The back half is colored dark blue to indicate 25 Gauge. Use with ET8198H StableGrip™ Vitreoretinal Handle (sold separately). Overall length 56 mm.



Strabismus

27041

Strabismus Scissors

Short blades. Straight. Blunt tips.
Overall length 110 mm.



27051

Strabismus Scissors

Long blades. Straight. Blunt tips.
Overall length 103 mm.



27052

Strabismus Scissors

Long blades. Curved.
Overall length 103 mm.



27070

Strabismus Scissors, Mod. Bonn

Curved. Overall length 88 mm.



E35701

Straight Strabismus Scissors

Blunt tips. 29 mm from mid-screw to tip.
Overall length 109 mm.



E35711

Curved Strabismus Scissors

Blunt tips. 29 mm from mid-screw to tip.
Overall length 109 mm.



Scissors

Strabismus

E35721

Curved Strabismus Scissors

Blunt tips. Strongly Curved. 27 mm from mid-screw to tip. Overall length 106 mm.



E3577

Knapp Curved Strabismus Scissors

Blunt tips. Lightly Curved. 25 mm from mid-screw to tip. Overall length 118 mm.



E3578

Knapp Curved Strabismus Scissors

Blunt tips. Strongly Curved. 25 mm from mid-screw to tip. Overall length 118 mm.



Suture Stich

E3321

Westcott Type Stich Scissors

Sharp pointed tips. 21 mm from mid-screw to tip. Overall length 120 mm.



E3340

STORZ Stich Scissors Fine Straight

Sharp pointed tips. 17 mm from mid-screw to tip. Flat serrated handle with polished finish. Overall length 112 mm.



Tenotomy

27012

Stevens Tenotomy Scissors

Short. Curved. Overall length 103 mm.



27060

Westcott Tenotomy Scissors

Right. Blunt/blunt. Slightly curved.
Overall length 110 mm.



27062

Westcott Tenotomy Scissors

Right. Sharp pointed/sharp pointed.
Overall length 110 mm.



27082

Stevens Tenotomy Scissors


Long. Curved. Overall length 110 mm.



27100

Westcott Micro Tenotomy Scissors

Blunt/Blunt. Overall length 99 mm.



E3320L

Wetscott Curved Tenotomy Scissors

Left. Blunt tips. 19 mm from mid-screw to tip.
Overall length 115 mm.



E3320R

Wetscott Curved Tenotomy Scissors

Right. Blunt tips. 19 mm from mid-screw to tip.
Overall length 115 mm.



Scissors

Tenotomy

E3320RS

Shepard-Westcott Tenotomy Scissors

Right. Blunt tips. Curved. 21 mm from mid-screw to tip.
Overall length 123 mm.



E3320RX

Westcott Curved Tenotomy Scissors

Right. Sub-mini model. Blunt tips. 9 mm from mid-screw to tip.
Overall length 88 mm.



E3552I

Stevens Curved Tenotomy Scissors

Blunt tips. 19 mm from mid-screw to tip.
Overall length 98 mm.



E3554I

Stevens Straight Tenotomy Scissors

Light model. Blunt tips. 19 mm from mid-screw to tip.
Overall length 98 mm.



E3556I

Stevens Curved Tenotomy Scissors

Light model. Blunt tips. 19 mm from mid-screw to tip.
Overall length 98 mm.



E3560

Stevens Straight Tenotomy Scissors

Blunt tips. Overall length 112 mm.



E3562

Stevens Curved Tenotomy Scissors

Blunt tips. Overall length 112 mm.

